

extrostyle Medical Shoes



EXTROSTYLE

EXTROSTYLE IS A COMPANY IN CONSTANT EVOLUTION, ANCHORED TO ITS PHILOSOPHY BUT ORIENTED TOWARDS CONTINUOUS DEVELOPMENT.

Wellness: the value that inspires all our actions, from the design to the creation of footwear. Craftsmanship and passion are the elements that allow us to make anti-fatigue shoes that increase reactivity in situations of risk; because we are firmly convinced that "WALKING IS LIVING".

Our footwear reduces fatigue and improves the quality of life of those who choose to face life with a smile, aware of their skills and their value. That is why for over thirty years we have accompanied our customers in daily life, including work, making every movement pleasant. In our factories, located in the **Umbria** and the **Marche** region, attention to detail and a **100% Made in Italy** production chain are combined with an orientation towards **innovation in the footwear market** and **environmental sustainability**, our strengths.

Every day we demand the highest quality in every production process: all ExtroStyle footwear has the CE mark for **Medical Devices** (MDR directive 745/17) as class I surgical medical devices registered at the Ministry of Health, as well as the UNI EN ISO 9001 and UNI EN ISO 13485 **system certifications**.

OUR DISTRIBUTORS







CERTIFICATIONS AND RELIABILITY

















All Extrostyle products have the **CE mark** in accordance with the MDR directive 745/17. To these, Extrostyle adds the UNI EN ISO 9001 and UNI EN ISO 13485 **System certifications**. The integrated management system allows to organise and monitor the company and its products at 360°.

Our organisation today allows us to offer a unique service, ensuring our customers delivery throughout Italy **WITHIN 24/48 HOURS** after receiving the order.

INDEX

INSOLE SYSTEM	5
ACTIVE LINE	7
ENERGY LINE	21
ENERGY 10	33
COMFORT LINE	35
REHABILITATION LINE	57
DIABETIC FOOT BRACES AND POST-OPERATIVE DEVICES LINE	67
EXTREME LINE	78





ORTHOSES PRESCRIPTION DEPENDING ON THE RISK PHASE

LEVEL OF DIABETIC RISK	PROBLEM	GUIDELINE FOR THE (FOOTWEAR AND OR		IDEAL PRODUCT
O LOW RISK	Diabetic patient with normal conformed foot. No neuropathy.	Calibrated FOOTWEAR Comfort	ORTHOPEDIC INSOLE No one or protective	Active Line Energy Line Comfort Line Rehabilitation Line
1 MEDIUM RISK	Neuropathy with light deformities.	Perfect FOOTWEAR Theraupedic to avoid stress caused by an inappropriate upper OUTSOLE Flexible bio-mechanic	ORTHOPEDIC INSOLE Protective and customized	Active Line Energy Line Rehabilitation Line
2 HIGH RISK	Neuropathy or arterial desease. Relevant deformities but without ulceration and/or previous amputation.	Accomodate d rebalance pr FOOTWEAR Theraupedic OUTSOLE Flexible bio-mechanic or semi-rigid	eformities and essure peaks ORTHOPEDIC INSOLE Customized	Active Line Energy Line Rehabilitation Line
3 VERY HIGH RISK	Neuropathy or arterial desease. Deformities, ulceration and/or previous amputation.		comodate deformities pressure peaks ORTHOPEDIC INSOLE Customized	Active Line Energy Line Diabetic foot and post operative device





EXTROSTYLE INSOLE SYSTEM

The insoles used for the ACTIVE and ENERGY lines are composed of a high-density open-cell foam that ensures **absorption and deabsorption** capacity, keeping the foot ventilated, dry and odourless, as well as ensuring a pleasant sensation when walking. The **ergonomic design** of the footbed offers high comfort, minimising the effort and tensions in the plantar support, both during walking and in a static position, dampening the impact of pressure during walking.

position, dampening the impact of pressure during walking. CHARACTERISTICS OF THE ENERGY LINE

→3 mm footbed

INSOLES

To be used immediately after an operation on the forefoot in the presence of very large dressings and bandages.

→6 mm footbed

To be used during rehabilitation and in the presence of bandages and swellings of the forefoot.

→9 mm footbed

To be used as a standard footbed already inside the footwear.



arneplant[®]

FIBRE INSOLE SYSTEM

Fibre insoles and TECNICA footwear are great allies against diabetic foot issues, as they allow the authorised healthcare professional to use this footwear at all stages of risk, respecting all international protocols for the treatment and care of diabetic foot issues.

CHARACTERISTICS OF FIBRE INSOLES

S-Rigid insole

Extra-rigid fibreglass insole. It makes the outsole completely rigid, limiting the joint stress of the foot during walking and reducing significantly the pressure peaks on the forefoot, redistributing them correctly.

S-Semirigid insole

Semi-rigid fibreglass insole. It makes the outsole semi-rigid, supporting the foot without forcing it. The pressure peaks on the forefoot are reduced and the load on the joints is limited.







S-CLOUD INSOLE SYSTEM

The S-Cloud insole is supplied as an accessory for the ENERGY, REHABILITATION and Diabetic foot and post operative device lines, and for the Tecnica Stelt and Tecnica Biomed products. The innovative S-Cloud insole allows to remove little squares in the lower part of the entire insole surface thanks to the pre-incision of the 6 mm Evaform. The orthopaedic technician will identify the pieces to be removed to allow the

immediate discharge of pressures in the ulcerated area or in the area at risk of ulceration, allowing to create a precise zone of unloading of the necessary size in a precise point. The 3 mm self-modeling memory foam inside the insole, reduces the forces of the edges of the unloading area on the plantar area, distributing thus the body weight over a wider surface avoiding in this way pressure peaks.

CHARACTERISTICS OF S-CLOUD INSOLES


- **1**. Breathable, antibacterial, antimicrobial, antiodour and hypoallergenic On-Steam cover.
- 2. 3 mm self-moulding memory foam.
- 3. 10 mm x 10 mm modular inserts in Evaform, 6 mm
- **4.** Fibre adhesive base to stabilise the Evaform modules.

COMBINATION OF USE

The S-Cloud plantar system is divided into five lines of use:

- →S-Cloud Energy compatible with the ENERGY line (pair)
- →S-Cloud Fisio compatible with the REHABILITATION line (pair)
- →S-Cloud Diab

 compatible with the Diabetic foot braces line
 (piece)
- →S-Cloud Stelt (piece)
- →S-Cloud Biomed (piece)





S-AIR INSOLE SYSTEM

The S-Air insoles are included throughout the REHABILITATION LINE. They provide a pleasant

sensation when walking, avoiding also the consequences of the continuous impact on the joints.

CHARACTERISTICS OF S-AIR INSOLES

- → Full ventilation thanks to the 100% breathable foam, which allows the air circulation.
- → Keeps feet dry eliminating humidity through its absorption system.
- →Eliminates odours, neutralising germs thanks to an antibacterial treatment.
- → Prevents fatigue thanks to its shock absorption capacity.
- → Longer-lasting foam thanks to its storage capacity.







ACTIVE LINE

ACTIVE is a line of footwear designed specifically to protect sensitive feet during physical activity. The biomechanical ACTIVE outsole allows a natural transition between the three phases of walking, ensuring a reduction of the pressure peaks on the forefoot and limiting the wear on the joints, thus allowing a more spontaneous walk and reducing the feeling of fatigue. The self-moulding and seamless

upper at the rubbing points envelops the anatomy of the foot and its possible deformities, ensuring comfort and preventing painful friction. The footwear of the ACTIVE line combined with fibre insoles can cover all stages of risk as they follow international protocols for the treatment and care of diabetic foot issues.

DEFINITION OF THE OUTSOLE RIGIDITY

- → ACTIVE FLEXIBLE: is when inside the shoe there is placed only the Tecnica insole.
- → ACTIVE SEMI-RIGID: is when the semi-rigid fiber glass insert is placed inside the shoe + the Tecnica insole.
- →ACTIVE RIGID: is when the rigid fiber glass insert is placed inside the shoe + the Tecnica insole.

DESCRIPTION

- →ACTIVE FLEXIBLE: indicated for the treatment of rheumatic foot and for the prevention of grade 0 and 1 diabetic foot in the presence of neuropathies accompanied by vascular diseases. The flexible biomechanical outsole facilitates the footstep rolling and distributes the pressure peaks.
- → ACTIVE SEMI-RIGID: indicated for the treatment of rheumatic foot suffering from pain, and the prevention of grade 1 and 2 diabetic foot. The semi-rigid biomechanical sole is also ideal for diabetic patients who need maximum protection against ulcer risk.
- →ACTIVE RIGID: indicated for the prevention of third level diabetic foot with the need to relieve the plantar pressure, designed for those who suffer from previous ulcerations or amputations, also indicated for patients who suffer from rheumatism and have serious deformities. The rigid biomechanical outsole limits the joint stress of the foot during walking and significantly reduces the pressure peaks on the forefoot and redistributes them correctly.

PLANTARE S-RIGID ACTIVE

Certification no. 2080794

Sizes from 35 to 46

Ref. no. EXO51150

PLANTARE S-SEMIRIGID ACTIVE

Certification no. 2080795

Sizes from 35 to 46

Ref. no. EXO51150

PLANTARE ACTIVE 6 mm

Certification no. 2142381

Sizes from 35 to 46

Ref. no. EX050530









ACTIVE PREVENTION FOOTWEAR LINE

















The **PADDED LINING**, which is comfortable, breathable and **with no seams**, welcomes the foot in a protective and cushioning casing.

Tecnica insole, removable and in microfibre

PADDED LINING



The central part of the outsole, being wide, guarantees **real stability** with consequent relaxation for the entire musculature. In addition, the Active outsole allows a natural "rolling" of the foot, without brusque movements.















STITCHING








SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

LINING











EXTRA FLEXIBLE ANCHE SHOCK ABSORBING ACTIVE OUTSOLE



PADDED AND SOLID ACTIVE HEEL SUPPORT



















SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

INSOLE

AND SHOCK ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

STITCHING

DEVICE

LINING

ACTIVE ALBA BIANCO DONNA

Certification no. 2258605

Color	White
Sizes	from 35 to 43
Ref. no.	EX055400





ACTIVE SAVONA ARGENTO DONNA

Certification no. 2258606

Color	Silver
Sizes	from 35 to 43
Ref. no.	EX055400






ACTIVE GELA BIANCO/ROSSO DONNA

Color	White/Red
Sizes	from 35 to 43
Ref. no.	EX055200













EXTRA FLEXIBLE AND SHOCK ABSORBING OUTSOLE













PADDED LINING



















REMOVABLE EXTRAF INSOLE AND S

AND SHOCK
ABSORBING OUTSOLE

PREDISPOSITION NO INNER FOR INSOLE STITCHING

MEDICAL DEVICE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

Color

Sizes

LINING

Black/Orange

from 35 to 43

ACTIVE CATANIA NERO DONNA

Certification no. 2021487



ACTIVE CAPRI GRIGIO/NERO DONNA

Certification no. 2080428

Color	Grey/Black
Sizes	from 35 to 43
Ref. no.	EX055100



ACTIVE PALERMO BIANCO DONNA

Certification no. 2021485

Color	White/Blue
Sizes	from 35 to 43
Ref. no.	EX054900





PADDED AND SOLID ACTIVE HEEL SUPPORT









FOR INSOLE









DEVICE PREVENTION OF DIABETIC FOOT

PADDED LINING









TOTAL OPENING WITH DOUBLE ZIP



EXTRA FLEXIBLE ANCHE SHOCK ABSORBING ACTIVE OUTSOLE









AND SHOCK

ABSORBING OUTSOLE





STITCHING



DEVICE





LINING DIABETIC FOOT

ACTIVE MARATEA NERO DONNA

Certification no. 2020302





Color	Black
Sizes	from 35 to 43
Ref. no.	EX056100

ACTIVE ENNA NERO/ORO DONNA

Certification no. 2021488





Color	Black/Gold Glitter
Sizes	from 35 to 43
Ref. no.	EX055300

ACTIVE TROPEA NERO/ORO DONNA





Color	Black/Gold Glitter
Sizes	from 35 to 43
Ref. no.	EX055500











ABSORBING OUTSOLE











FOR INSOLE STITCHING

MEDICAL DEVICE PREVENTION OF DIABETIC FOOT

LINING

ACTIVE MILAZZO GRIGIO UOMO

Certification no. 2142374

Color	Black/Grey
Sizes	from 39 to 46
Ref. no.	EX055200

ACTIVE PALERMO BIANCO UOMO

Color White/Blue Sizes from 39 to 46 Certification no. 2021476 Ref. no. EX055200





ACTIVE PALERMO NERO UOMO

Certification no. 2258590

Color	Black/Bordeaux
Sizes	from 39 to 46
Ref. no.	EX055200

ACTIVE ELIANO BLU UOMO

Certification no. 2021482

Color	Blue/Grey
Sizes	from 39 to 46
Ref. no.	EX055500





ACTIVE STILO GRIGIO UOMO

Color	Grey
Sizes	from 39 to 46
Ref. no.	EX055600











AND SHOCK

ABSORBING OUTSOLE



FOR INSOLE



STITCHING







MEDICAL SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

PADDED LINING

ACTIVE ELIA NERO UOMO

Certification no. 2021481



Color Black/White
Sizes from 39 to 46
Ref. no. EX055500

ACTIVE CATANIA NERO UOMO

Certification no. 2021478



Color Black/Dark brown
Sizes from 39 to 46
Ref. no. EXO55300

ACTIVE TRAPANI 14 T.MORO UOMO



























REMOVABLE INSOLE

AND SHOCK ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

MEDICAL DEVICE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

PADDED LINING

ACTIVE ALCAMO NERO UOMO









TOTAL OPENING WITH DOUBLE ZIP



EXTRA FLEXIBLE ANCHE SHOCK ABSORBING ACTIVE OUTSOLE





















LINING

REMOVABLE **INSOLE**

AND SHOCK ABSORBING OUTSOLE

FOR INSOLE

STITCHING

MEDICAL **DEVICE**

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

ACTIVE TARANTO NERO UOMO

Certification no. 2121483

Color Black/ Bordeaux Sizes from 39 to 46



ACTIVE CONERO NERO UOMO

Certification no. 2142376


Color Black Sizes from 39 to 46 Ref. no. EX056300



INTERNAL MEASURES












 $\textbf{MEASURES OF THE MODELS:} \ \textbf{FORMIA}, \ \textbf{ALBA}, \ \textbf{SAVONA}, \ \textbf{GELA}, \ \textbf{VIESTE}, \ \textbf{CATANIA}, \ \textbf{CAPRI}, \ \textbf{PALERMO}, \ \textbf{MARATEA}, \ \textbf{SCALEA}, \ \textbf{ENNA}, \ \textbf{TROPEA}, \ \textbf{TR$ MILAZZO, ELIANO, STILO, ELIA, ALCAMO, TARANTO, CONERO.

				FI	T 12								
Sizes	UM	35	36	37	38	39	40	41	42	43	44	45	46
Lenght	mm	238	244	251	258	264	271	278	284	291	298	304	311
Forefoot width	mm	80	81,5	83	84	86	87,5	89	93	94	95,5	96,5	98
Heel width	mm	56	57	57,5	59	60	61	62	65	66	67	68	69
Forefoot circumference	mm	231	235,5	239,5	244	249,5	253	258	268	272	277	281	286
Toe height	mm	29	29,5	30	30,5	31	31,5	32	32,5	34	34,5	35	35,5
Instep circumference	mm	231	235,5	240	244,5	249	253,5	258	267,5	272	276,5	281	285,5

MEASURES OF THE MODELS: ACTIVE TRAPANI TESTA MORO LIOMO

FIT 14													
Sizes	UM	35	36	37	38	39	40	41	42	43	44	45	46
Lenght	mm	238	244	251	258	264	271	278	284	291	298	304	311
Forefoot width	mm	83	84	85,5	87	88,5	90	92	95	97	99	100	102
Heel width	mm	57	58	59	60	61	62	63	66	67	69	69	70
Forefoot circumference	mm	241	245	251	255	260	264	268	280	284	291	294	298
Toe height	mm	31,5	32	32,5	33	33,5	34	34,5	37	37,5	37,5	38	39
Instep circumference	mm	241	245,5	250	254,5	259	263,5	268	280	284,5	289	293,5	298





ENERGY LINE

From the technology and the experience of ExtroStyle, and from scientific studies and advanced researches in the biomechanical field, results the new ENERGY line, a performant footwear with a current and innovative design in many respects, suitable for

the treatment of foot morphology alterations. The footwear of the ENERGY line combined with our fibre reinforcements can cover all stages of diabetic foot risk as they follow the international protocols for the treatment and care of diabetic foot issues.

DEFINITION OF THE OUTSOLE RIGIDITY

- → ENERGY FLEXIBLE is when inside the shoe there is placed only the Tecnica insole.
- → ENERGY SEMI-RIGID is when the semi-rigid fiber glass insert is placed inside the shoe + the Tecnica insole.
- → ENERGY RIGID is when the rigid fiber glass insert is placed inside the shoe + the Tecnica insole.

DESCRIPTION

- → ENERGY FLEXIBLE: is indicated for the treatment of the rheumatic foot disease and for the prevention of diabetic foot with risk degrees 0 and 1 in presence of neuropathies accompanied by vasculopathies. The flexible biomechanical outsole facilitates the footstep rolling and distributes pressure peaks.
- → ENERGY SEMI-RIGID: is indicated for the treatment of rheumatic foot suffering from pain, and for the prevention of diabetic foot with risk degree 1 and 2. The biomechanical semi-rigid outsole is also ideal for diabetic patients who need maximum protection against the risk of ulceration.
- → ENERGY RIGID: is indicated for the third level of diabetic foot prevention, with the need to relieve plantar pressure, designed for those suffering from previous ulcerations or amputations, also suitable for patients suffering from rheumatism or patients with severe foot deformities. The rigid biomechanical outsole limits the foot joint stress during walking and reduces significantly pressure peaks on the forefoot and redistributes them correctly.

PLANTARE ENERGY

Certification no. 3 mm 1766123 6 mm 1766124 9 mm 1766125 Sizes from 35 to 50

Ref. no. EX030650

PLANTARE S-RIGID ENERGY

Certification no. 1754064

Sizes from 35 to 50

Ref. no. EXO31150





PLANTARE S-CLOUD ENERGY

Certification no. 2142386

Sizes from 35 to 50

Ref. no. EX031490

PLANTARE S-SEMIRIGID ENERGY

Certification no. 2142379

Sizes from 35 to 50

Ref. no. EXO31150









ENERGY LINE VITAL ENERGY AT YOUR FEET

















FOR THE TREATMENT OF ALTERATIONS IN THE MORPHOLOGY OF THE FEET



















REMOVABLE INSOLE

EXTRA FLEXIBLE AND SHOCK ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

MEDICAL DEVICE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

PADDED LINING

ENERGY PERLA

Certification no. 2421895

Color	Grey/ Cream
Sizes	from 35 to 43
Ref. no.	EX035600





REMOVABLE TECNICA INSOLE

ENERGY AGATA

Certification no. 2421896

Color	Grey/ Cream
Sizes	from 35 to 43
Pof no	EV03E/00







ADJUSTABLE VELCRO OPENING

ENERGY PETRA

Color	Dove grey
Sizes	from 35 to 43
Ref. no.	EX035200





BIOMECHANICAL OUTSOLE



















INSOLE

EXTRA FLEXIBLE AND SHOCK ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

MEDICAL DEVICE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

LINING

ENERGY MARINA

Certification no. 1903605

Color Black
Sizes from 35 to 43

Ref. no. EX035100





ADJUSTABLE DOUBLE VELCRO OPENING

ENERGY ALICE

Certification no. 1717176

Color Bronze
Sizes from 35 to 43
Ref. no. EX034900





COMFORTABLE, HARMONIOUS AND ELEGANT



















LINING

AND SHOCK ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

MEDICAL DEVICE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

ENERGY VALENTINA

Certification no. 1903609





Color	Black
Sizes	from 35 to 43
Ref. no.	EX035700

ENERGY ELENA


Certification no. 1765781





Color	Black
Sizes	from 35 to 43
Ref. no.	EX034600

ENERGY LUCIA





Color	Black
Sizes	from 35 to 43
Ref. no.	EX035200





Color Grey/Black Sizes from 39 to 50

Ref. no. EX035600







REMOVABLE INSOLE



AND SHOCK ABSORBING OUTSOLE



PREDISPOSITION FOR INSOLE



NO INNER STITCHING



MEDICAL DEVICE



SUITABLE FOR THE PREVENTION OF DIABETIC FOOT



PADDED LINING







REMOVABLE INSOLE



EXTRA FLEXIBLE PI AND SHOCK ABSORBING OUTSOLE



PREDISPOSITION FOR INSOLE



NO INNER STITCHING



MEDICAL DEVICE



SUITABLE FOR THE PADDED PREVENTION OF LINING DIABETIC FOOT



ENERGY LORIS

Certification no. 2421682

Color	Black
Sizes	from 39 to 50
Ref. no.	EX035300





BREATHABLE UPPER

ENERGY LUIGI

Certification no. 1903607

Color	Black
Sizes	from 39 to 50
Ref. no.	EX035300





BIOMECHANICAL OUTSOLE

ENERGY VALENTINO NERO

Color	Black
Sizes	from 39 to 50
Ref. no.	EX035600





BIOMECHANICAL OUTSOLE







INSOLE

EXTRA FLEXIBLE AND SHOCK ABSORBING OUTSOLE



FOR INSOLE

NO INNER

STITCHING



MEDICAL SUITABLE FOR THE PREVENTION OF DIABETIC FOOT



LINING

ENERGY PAOLO MARRONE

MARRONE	
Certification no.	2060235 _

Color	Brown
Sizes	from 39 to 50
Ref. no.	EX035300

ENERGY PAOLO NERO

Certification no. 2060235

Color	Black
Sizes	from 39 to 50
Ref. no.	EX035300





ENERGY CARLO

Certification no. 1765782

Black
from 39 to 50

Ref. no. EX035600

ENERGY DANIELE

Certification no. 1696618

Color	Black	
Sizes	from 39 to 50	






ENERGY ALEX NERO

Color	Black
Sizes	from 39 to 50
Ref. no.	EX035100















PHASE 1 PREVENTION

Thanks to the predisposition for the Energy 10 **removable insole**, it is a great ally for the prevention or containment of the most common foot diseases such as: hallux valgus, hammertoes, metatarsalgia, heelitis and vascular diseases as well as for the prevention with a 0 - 1 risk in diabetic foot. The **flexible biomechanical energy outsole** allows a natural transition between the three phases of walking.

ensuring a reduction of the pressure peaks on the forefoot and limiting the wear on the joints, thus allowing a more spontaneous walk and reducing the feeling of fatigue. The **self-moulding and seamless upper** at the rubbing points envelops the anatomy of the foot and its possible deformities, ensuring comfort and preventing painful friction.













HALLUX VALGUS HAMMER TOES
CLAW TOES
METATARSALGIA

BURSITIS CHARCOT TARSAL TUNNEL

PHASE 2 POST-OPERATIVE REHABILITATION

The energy 10 **full opening shoe with 2 zips**, equipped with a flexible outsole, can be defined as a multifunctional shoe thanks to a combination of insoles that are variable in height. **It can be worn from the first day after the operation until the end of convalescence**. This implies a more correct and balanced walking and distribution of body weight, as well as an undeniable benefit for the back, pelvis and legs. Also the psychological well-being is not to be underestimated, since the energy 10 footwear can be worn on both feet


and immediately after the operation. A **removable heel** (patented system) inclined by 8° can be applied to keep the foot in talism for all rehabilitation courses where the surgeon deems it necessary. The Energy 10 footwear was subjected by the Moticon Science (independent company) of Munich, to test the **PLANTAR DISCHARGE**, comparing it with five other footwear of other brands among the most sold on the European market and was the **BEST AND MOST EFFECTIVE IN ALL THE COMPARISON TESTS**.

TOTAL OPENING THANKS TO THE 2 ZIPS DIFFERENT HEIGHTS OF THE INSOLES ACCORDING TO THE NEEDS

REMOVABLE HEEL PATENTED TO KEEP
THE FOOT IN TALISM







ENERGY 10, THE BEST EUROPEAN FOOTWEAR FOR PLANTAR DISCHARGE!

Study conducted by Moticon Science of Munich

PHASE 3 PROTECTION OF FOOT HEALTH

The Energy 10 footwear has been designed for men and women who want to live their natural "vital energy". The innovative technical characteristics open the way to a new idea of well-being, health and vitality. Indeed the line

is designed to meet the needs and to prevent and manage the problems of the sensitive foot but above all, to live well.







REMOVABLE EXTRA FLEXIBLE **INSOLE** AND SHOCK ABSORBING OUTSOLE



PREDISPOSITION FOR INSOLE





NO INNER STITCHING



MEDICAL **DEVICE**



SUITABLE FOR THE PREVENTION OF DIABETIC FOOT



PADDED LINING



ENERGY 10 BLU

Certification no. 2421679

Color	Blue
Sizes	from 35 to 41 (woman)
	from 42 to 50 (man)

Ref. no. EX034750

ENERGY 10 JEANS

Certification no. 2060235

Color Jeans/Blue Sizes from 35 to 41 (woman) Ref. no. EX035070





ENERGY 10 BIANCO

Certification no. 2056718

Color	White / Black							
Sizes	from 35 to 41 (woman)							
Ref. no.	EX034750							



Е

С

NERGY PLUS	Color	Black
Certification no. 1766122	Sizes	from 35 to 46
	Ref. no.	EX031050
	Size TABLE	
	35-36	
E Common of the	37-38	
The same of the sa	39-40	
	41-42-43	
AL.	44-45-46	
######################################		100





UNISEX









FXTRA FI FXIBI F AND SHOCK ABSORBING OUTSOLE



PREDISPOSITION FOR INSOLE



NO INNER STITCHING



MEDICAL **DFVICE**



PADDED LINING

PREVENTION OF DIABETIC FOOT

ENERGY MORRIS RIGID

Certification no. 1903608

Color Blue from 35 to 50 Sizes Ref. no.

EX035600

ENERGY JACQUES FLEX

Certification no. 1920242

Color Blue Sizes from 35 to 50 Ref. no. EX033900






The ENERGY MORRIS RIGID and ENERGY JACQUES FLEX, thanks to the double velcro closure, have a complete opening in the anterior part of the shoe, making it more practical to wear. The extra wide fit and the adjustable internal volume make these shoes suitable for patients with large, bandaged, swollen or edematous feet. In addition, they are designed for the insertion of a customized insole.

→ SEMI-RIGID OUTSOLE: to be used in primary

prevention with risk degree 0 or 1, polyarthritis, Morton's neuroma, osteoarthritis, etc.

→ RIGID OUTSOLE: to be used in secondary prevention with risk degree 2 and 3, for patients with previous ulcerations at risk of recurrence, for patients with deformities, ulcerations or previous amputations, arthrodesis of the toes, arthrodesis of the tarsus, Charcot foot or biomechanical problems of movement.

INTERNAL MEASURES









MEASURES OF THE MODELS: PERLA, AGATA, PETRA, VALENTINA, MARINA, ELENA, LUCIA, ALICE, IVAN, VALENTINO, LORIS, LUIGI, PAOLO, CARLO, DANIELE, ALEX, ENERGY 10, MORRIS, JACQUES.

Sizes	UM	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50
Lenght	mm	233	240	247	253	260	267	273	280	287	293	300	307	313	320	327	333
Forefoot width	mm	84,5	86	87,5	89	90,5	92	93,5	96,5	98	99,5	101	102,5	104	105,5	107	108,5
Heel width	mm	56,4	57,6	58,8	60	61,5	62,4	63,6	66,8	68	69,2	70,4	71,6	72,8	74	75,2	76,4

INSOLE 3 mm																	
Forefoot circumference	mm	232,5	237	241,5	246,5	251	255,5	260	267	271,5	276	280,5	285	289,5	294	298,5	302,5
Toe height	mm	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42

INSOLE 6 mm																	
Forefoot circumference	mm	227	231,5	236	240,5	245	249,5	254	261	265,5	270	274,5	279	283,5	288	292,5	297
Toe height	mm	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39

INSOLE 9 mm																	
Forefoot circumference	mm	221	225,5	230	234,5	239	243,5	248	255	259,5	264	268,5	273	277,5	282	286,5	291
Toe height	mm	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36



COMFORT LINE

COMFORTABLE, SOFT, SOLID AND RESISTANT FOOTWEAR

The COMFORT line consists of a collection of unique, functional and practical footwear which combine technical comfort and style in the same product. They are suitable for all ages and all seasons. All models have a large internal volume for the foot support and are suitable for the insertion of custommade insoles. The trendy colours and the lining in On Steam® microfibre material offer a very pleasant

touch sensation, guaranteeing 100% breathability. The heel, of ideal height, is designed to guarantee a correct posture and offer a right boost for walking. The attention to detail, the modern design and the constant research for materials make TECNICA a product of great quality, guaranteed by an Italian company that for over thirty years has designed and manufactured products of excellence.





ON STEAM®

THE WHOLE COMFORT LINE HAS THE INSOLE AND THE SHOE LINING COVERED BY ON STEAM® MICROFIBER, WHICH ABSORBS ODORS, IS ANTIBACTERIAL AND BREATHABLE.



On Steam® is the new 100% breathable microfiber, with a sensation like leather and with a capacity of humidity absorption 8 times its weight in water. It provides the maximum comfort, total dryness, and a thermo-regulating effect. Sweat can not accumulate inside the footwear because the textile leads it out, where it evaporates, maintaining thus the foot dry in every moment. The main characteristics of this microfiber are: absorption action, desorption action,

maximum water removal, complete dryness, high abrasion resistance, it is antibacterial, antimicrobial, and odorless, with a thermo-regulating effect, it is also hypoallergenic, thermoformable, and Chromium VI-free, and has an ecological production system without CO2. Furthermore, this microfiber has also the OEKO -TEX Standard 100 ecological certification: an external laboratory has submitted the textile to strict tests and has certified the total absence of harmful health agents.

PLANTARE COMFORT

Certification no. 2473528

Sizes from 35 to 43

Ref. no. EX040650

PLANTARE CORK MEMORY

Certification no. 2473529

Sizes from 35 to 42

Ref. no. EX040700





PLANTARE CORK SLOW

Certification no. 2473530

Sizes from 39 to 50

Ref. no. EX040750



THE PLANTARE **COMFORT** CAN BE INSERTED IN: TECNICA ADELE PETROLIO, TECNICA ADELE ROSSO, TECNICA ADELE OCRA, TECNICA ADELE SUDAN, TECNICA ADELE FIORI, TECNICA ANNA FIORI, TECNICA SIMONA LILLA, TECNICA SIMONA VERDE, TECNICA SIMONA FIORI, TECNICA SIMONA MERLOT, TECNICA SIMONA BIANCO/PLATINO, TECNICA SIMONA GRIGIO/SABBIA, TECNICA SONIA ROSSO, TECNICA SONIA BLU, TECNICA GRAZIA CELESTE, TECNICA GRAZIA SABBIA, TECNICA GRAZIA NERO/GRIGIO, TECNICA CLAUDIA LAMINATO ROSSO, TECNICA KIARA BLU, TECNICA KIARA GHIACCIO PERLA/ARGENTO, TECNICA BEA GRIGIO, TECNICA ALMA BLU, TECNICA CARLA, TECNICA GEA, TECNICA BETTA, TECNICA GAIA LANA ARANCIO, TECNICA GAIA LANA ROSSO, TECNICA GAIA LANA BLU, TECNICA GAIA LANA BLU, TECNICA GAIA LANA BLU, TECNICA GAIA LANA PERDE, TECNICA ALEXIA.

THE PLANTARE **CORK MEMORY** CAN BE INSERTED IN: TECNICA CRISTINA ARGENTO, TECNICA ANGELA SABBIA, TECNICA ANGELA NERO, TECNICA MARY BORDEAUX, TECNICA MARY VERDE, TECNICA SARA LILLA, TECNICA SARA FIORI, TECNICA MARY NERO, TECNICA MARY SABBIA, TECNICA MARY BRONZO, TECNICA TANIA, TECNICA LOREN, TECNICA DOROTY NERO, TECNICA DORIS GRIGIO, TECNICA ASIA, TECNICA BENEXA.

THE PLANTARE CORK SLOW CAN BE INSERTED IN: TECNICA RICCARDO NERO, TECNICA RICCARDO, TECNICA RICK NERO, TECNICA RICK, TECNICA BENEXA, TECNICA ANGELO, TECNICA DEVIS, TECNICA DYLAN.























PADDED LINING







INSOLE











EXTRA FLEXIBLE PREDISPOSITION
AND SHOCK FOR INSOLE
ABSORBING OUTSOLE

NO INNER STITCHING

MEDICAL DEVICE

PADDED LINING

HORIZON CAORLE

Certification no. 2469202

Color	Brown
Sizes	from 35 to 43
Ref. no.	EX045700





ADJUSTABLE ELASTIC AND VELCRO OPENING

HORIZON ESTE

Certification no. 2469203

Color	Brown
Sizes	from 35 to 43
Ref. no.	EX045800





TOTAL OPENING

HORIZON JESOLO

Color	Black
Sizes	from 35 to 43
Ref. no.	EX045700





ZIP OPENING







REMOVABLE

INSOLE





PREDISPOSITION FOR INSOLE



NO INNER STITCHING



PADDED LINING



Certification no. 2469205

Color	Sand
Sizes	from 35 to 43
Ref. no.	EX045700





ADJUSTABLE ELASTIC AND VELCRO OPENING

HORIZON CADORE

Certification no. 2469206

Color	Brown
Sizes	from 35 to 43
Ref no	EX0/5800





ELASTIC UPPER NEAR THE BIG TOE

HORIZON GARDA

Color	Beige
Sizes	from 35 to 43
Ref. no.	EX045700





ELASTIC UPPER NEAR THE BIG TOE











AND SHOCK ABSORBING OUTSOLE FOR INSOLE



STITCHING



MEDICAL DEVICE PADDED LINING




INTERNAL MEASURES





Lenght







 $\label{lem:measures of the model:} \mbox{HORIZON CANAZEI - HORIZON CAORLE - HORIZON ESTE - HORIZON JESOLO - HORIZON CADORE - HORIZON GARDA.}$

Sizes	UM	35	36	37	38	39	40	41	42	43
Lenght	mm	234,6	241,3	248	254,6	261,3	267,9	274,6	281,3	287,9
Forefoot circumference	mm	223	227	232	236	241	245	250	254	259
Instep circumference	mm	225	229	234	238	243	247	252	256	261
Forefoot width	mm	81	82,5	84	85,5	87	88,5	90	91,5	93
Hell width	mm	52	53	54	55	56	57	58	59	60
Toe height	mm	18	19	20	21	22	23	24	25	26





















AND SHOCK FO
ABSORBING OUTSOLE

NO INNER STITCHING

MEDICAL DEVICE

PADDED LINING













STITCHING



DEVICE

PADDED LINING

TECNICA ADELE ROSSO

Certification no. 2421181

Color	Red
Sizes	from 35 to 43
Ref. no.	EX042900



TECNICA ADELE OCRA

Certification no. 2258636

Color	Ocher
Sizes	from 35 to 43
Ref. no.	EX042900



TECNICA ADELE SUDAN

Color	Sudan
Sizes	from 35 to 43
Ref. no.	EX042900



















REMOVABLE INSOLE

E AND SHOCK
ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

MEDICAL DEVICE

PADDED LINING

TECNICA ADELE FIORI

Certification no. 2421183

Color Flowers
Sizes from 35 to 43
Ref. no. EX042900





ADJUSTABLE DOUBLE VELCRO OPENING

TECNICA ANNA FIORI

Certification no. 2421184

Color Flowers
Sizes from 35 to 43
Ref. no. EX042900





STRETCH SIDE INSERT FOR DIFFERENT PATHOLOGIES























PREDISPOSITION FOR INSOLE REMOVABLE INSOLE EXTRA FLEXIBLE AND SHOCK ABSORBING OUTSOLE

NO INNER STITCHING

PADDED LINING

















REMOVABLE INSOLE

AND SHOCK ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

MEDICAL DEVICE

LINING

TECNICA SIMONA VERDE

Color	Green
Sizes	from 35 to 43
Ref. no.	EX042540

TECNICA SIMONA FIORI

Certification no. 2421179

Color	Flowers
Sizes	from 35 to 43
Ref. no.	EX042540





TECNICA SIMONA MERLOT

Certification no. 2421180

Color	Red
Sizes	from 35 to 43
Ref. no.	EX042540

TECNICA SIMONA BIANCO/PLATINO

Certification no. 2012533

Color	White/Platinum
Sizes	from 35 to 43
Ref. no.	EX042540





TECNICA SIMONA GRIGIO/SABBIA

Color	Grey/Sand
Sizes	from 35 to 43
Ref. no.	FX042540









REMOVABLE EXTRA FLEXIBLE AND SHOCK



EXTRA FLEXIBLE PREDISPOSITION AND SHOCK FOR INSOLE ABSORBING OUTSOLE



NO INNER STITCHING



MEDICAL DEVICE



LINING

TECNICA SONIA ROSSO

Certification no. 2258640

Color	Red
Sizes	from 35 to 43
Ref. no.	EX042540




Certification no. 2258641

Color	Blue
Sizes	from 35 to 43
Ref. no.	EX042540





Certification no. 2421192





Color Light blue/ Sand
Sizes from 35 to 43

Ref. no. EX042600

TECNICA GRAZIA SABBIA

Certification no. 2012532

Color	Dove grey/
	Sand
Sizes	from 35 to 43
Ref. no.	EX042600



TECNICA GRAZIA NERO/GRIGIO

Color	Black/Grey
Sizes	from 35 to 43
Ref no	EX0/3600



















REMOVABLE INSOLE

AND SHOCK
ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

R MEDICAL DEVICE

AL PADDED E LINING

TECNICA CLAUDIA LAMINATO ROSSO

Certification no. 24211930

Color	Red
Sizes	from 35 to 43
Ref. no.	EX042850



TECNICA KIARA BLU

Certification no. 1788338



Color	Blue
Sizes	from 35 to 43
Ref. no.	EX042750

TECNICA KIARA GHIACCIO PERLA/ARGENTO

Certification no. 2012535



Color Ice/Pearl/Silver
Sizes from 35 to 43
Ref. no. EX042750











AND SHOCK

FOR INSOLE ABSORBING OUTSOLE



NO INNER STITCHING



MEDICAL

DEVICE

PADDED LINING



Certification no. 2258643

Color Grey Sizes from 35 to 43 Ref. no. EX042750







SIDE INSERT



Color

Sand



TECNICA MARY BORDEAUX

Certification no. 2421119

Color Bordeaux Sizes from 35 to 42 Ref. no. EX042800





















REMOVABLE INSOLE

AND SHOCK
ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

STITCHING

MEDICAL DEVICE

PADDED LINING

TECNICA SARA LILLA

Certification no. 2421188



Color Lilac
Sizes from 35 to 42
Ref. no. EXO43060



STRETCH SIDE INSERT FOR DIFFERENT PATHOLOGIES

TECNICA SARA FIORI

Certification no. 2421189







STRETCH SIDE INSERT FOR DIFFERENT PATHOLOGIES







REMOVABLE EXTRA FLEXIBLE
INSOLE AND SHOCK
ABSORBING OUTSOLE



PREDISPOSITION FOR INSOLE



NO INNER STITCHING



MEDICAL

DEVICE

PADDED LINING



Certification no. 2056721

Color	Black
Sizes	from 35 to 42
Ref. no.	EX042800



TECNICA MARY SABBIA

Certification no. 2056720

Color	Sand
Sizes	from 35 to 42
Ref. no.	EX042800



TECNICA MARY BRONZO

Certification no. 1766119





Color Bronze Sizes from 35 to 42 Ref. no. EX042800

TECNICA TANIA























REMOVABLE INSOLE

EXTRA FLEXIBLE AND SHOCK ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

MEDICAL DEVICE

PADDED LINING









FOR INSOLE







EXTRA FLEXIBLE AND SHOCK ABSORBING OUTSOLE

STITCHING

DEVICE

PADDED LINING

TECNICA CARLA

Certification no. 2421185

Color	Grey
Sizes	from 35 to 43
Ref. no.	EX042950





ADJUSTABLE VELCRO OPENING

TECNICA GEA

Certification no. 2421186

Color	Beige
Sizes	from 35 to 43
Ref. no.	EX042950





TECNICA BETTA

Color	Blue
Sizes	from 35 to 43
Ref. no.	EX042950











REMOVABLE EXTRA FLEXIBLE INSOLE AND SHOCK ABSORBING OUTSOLE



PREDISPOSITION FOR INSOLE



NO INNER STITCHING



DEVICE

PADDED LINING



Certification no. 2045411

Color	Orange
Sizes	from 35 to 43
Ref. no.	EX042450





TECNICA GAIA LANA ROSSO

Certification no. 2045409

Color	Red
Sizes	from 35 to 43





Certification no. 2045408

Color	Blue
Sizes	from 35 to 43
Ref. no.	EX042450





TECNICA GAIA LANA VERDE

Color	Green
Sizes	from 35 to 43
Ref. no.	EX042450









REMOVABLE EXTRA FLEXIBLE INSOLE AND SHOCK



EXTRA FLEXIBLE PREDISPOSITION FOR INSOLE ABSORBING OUTSOLE



NO INNER STITCHING



MEDICAL

DEVICE

PADDED LINING





Color Blue
Sizes from 35 to 42
Ref. no. EXO42700



TECNICA DOROTY NERO

Certification no. 2203542

Color	Black
Sizes	from 35 to 42



TECNICA DORIS GRIGIO

Certification no. 2203540

Color	Grey
Sizes	from 35 to 42
Pof no	EV0/2000



TECNICA ASIA

Color	Black
Sizes	from 35 to 42
Ref. no.	EX042800









EMOVABLE EXTRA FLEXIBLE
INSOLE AND SHOCK
ABSORBING OUTSOLE



PREDISPOSITION FOR INSOLE



NO INNER MEDICAL STITCHING DEVICE



PADDED LINING

TECNICA RICCARDO NERO

Color	Black
Sizes	from 39 to 50
Ref. no.	EX043200

TECNICA RICCARDO

Certification no. 2123552

Color	Brown
Sizes	from 39 to 50
Ref. no.	EX043200





TECNICA RICK NERO

Certification no. 2418177

Color	Black
Sizes	from 39 to 50
Ref. no.	EX043500

TECNICA RICK

Certification no. 1754076

Color	Brown
Sizes	from 39 to 50
Ref no	EX0/3500





TECNICA BENEXA

Certification no. 1766118

Color	Black
Sizes	from 35 to 50
Ref. no.	EX043700



TECNICA ANGELO

Color	Blue
Sizes	from 39 to 46
Pof no	EV0/2000



















REMOVABLE INSOLE

EXTRA FLEXIBLE AND SHOCK ABSORBING OUTSOLE

PREDISPOSITION FOR INSOLE

NO INNER STITCHING

MEDICAL DEVICE

PADDED LINING

TECNICA DEVIS

Certification no. 2421194



 Color
 Blue

 Sizes
 from 39 to 50

 Ref. no.
 EX042600



ADJUSTABLE VELCRO CLOSURE

TECNICA DYLAN



Color	Beige
Sizes	from 39 to 50
Ref. no.	EX042600



REAL WOOL INSIDE AND ADJUSTABLE VELCRO CLOSURE







REHABILITATION LINE

A LINE DESIGNED FOR THE TREATMENT OF MOST FOOT DISEASES THAT COMBINES THE COMFORT AND SOFTNESS OF A SLIPPER WITH THE STABILITY AND RESISTANCE OF A SHOE.

This line of footwear is mainly used to improve walking in the elderly and in people who suffer from rheumatism. They are also used as post-intervention and rehabilitation footwear. The non-slip wide outsole is made with shock-absorbing compound that, combined with the thermoformable stabilising heel counters, guarantees a safer walk. The breathable upper is made of materials that shape around the forefoot and that allow to accomodate deformity or bandages without pain.

The velcro closure and the full opening upper allow those with joint mobility problems to comfortably wear the shoe with extreme ease, adjusting the velcro closure in the most optimal way. In addition, these footwear can be used in primary prevention with risk degree 0-1 in the diabetic foot, while by using the S-Rigid Fisio fibre insoles, which stiffen the outsole, it can be used in secondary prevention with risk degree 2 - 3.



































REMOVABLE EXTRA FLEXIBLE PREDISPOSITION INSOLE AND SHOCK FOR INSOLE ABSORBING OUTSOLE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

PADDED LINING

TECNICA MARISA

Certification no. 2470727



Color Cream Sizes from 35 to 42 Ref. no. EX043200

TECNICA MARTINA

Certification no. 2470728

Color	Cream
Sizes	from 35 to 42
Ref. no.	EX043200





Black

from 35 to 42

EX043200

TECNICA MIRIAM





















REMOVABLE EXTRA FLEXIBLE PREDISPOSITION INSOLE AND SHOCK FOR INSOLE ABSORBING OUTSOLE


STITCHING

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

PADDED LINING



Color Night blue Sizes from 35 to 50Ref. no. EX043300



TECNICA ANDRIA

























REMOVABLE EXTRA FLEXIBLE PREDISPOSITION INSOLE AND SHOCK FOR INSOLE ABSORBING OUTSOLE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

TECNICA AVOLA

Color	Black
Sizes	from 35 to 50
Ref. no.	EX044100



Certification no. 1919743

Color	Black
Sizes	from 35 to 50
Ref no	EX01/100





TECNICA SICILIA

Certification no. 2022819

Black
from 35 to 50

TECNICA DYNAMIC UOMO Certification no. 1787586

Color	Brown
Sizes	from 39 to 50
Ref. no.	EX013100





TECNICA GOLD 1

Certification no. 983653

Color	Black
Sizes	from 35 to 50
Ref. no.	EX013000



TECNICA GOLD AIR

Color	Grey
Sizes	from 35 to 50
Pof no	EV012000























REMOVABLE EXTRA FLEXIBLE PREDISPOSITION INSOLE AND SHOCK FOR INSOLE ABSORBING OUTSOLE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

TECNICA GOLD 3

Certification no. 983656

Color	Black
Sizes	from 35 to 50

Ref. no. EX013000



TECNICA FLEX

Certification no. 1304190

Color	Grey/Black
Sizes	from 35 to 50

Ref. no. EX012650

TECNICA DYNAMIC AIR

Certification no. 1765777

Color	Black
Sizes	from 35 to 50
Ref. no.	EX012370





TECNICA SOFT

Color	Black
Sizes	from 35 to 50
Ref. no.	EX012370























REMOVABLE EXTRA FLEXIBLE PREDISPOSITION INSOLE AND SHOCK FOR INSOLE ABSORBING OUTSOLE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

TECNICA 3T



Color Grey Sizes from 35 to 50 Ref. no. EX012260



VELCRO CLOSURE AND TOTAL OPENING

TECNICA 3S

Certification no. 1254501



Color Black Sizes from 35 to 50 Ref. no. EX012370





















PREVENTION OF DIABETIC FOOT

TECNICA 3M

Color	Blue
Sizes	from 35 to 50
Ref. no.	EX012260

TECNICA 3 AIR

Certification no. 1760034

Color	Grey
Sizes	from 35 to 50
Ref. no.	EX012430





TECNICA 3E

Certification no. 983660

Black
from 35 to 50

Ref. no. EX012370

TECNICA LARRY Certification no. 2142378

Color	Black
Sizes	from 35 to 50
Ref. no.	EX012430





TECNICA 3A

Certification no. 983657

Color	Grey
Sizes	from 35 to 50
Ref. no.	EX012430



TECNICA 1

Certification no. 983654

Certification no. 703034	Sizes	from 35 to 50
tecnica	Ref. no.	EX012700

Color

Grey























PADDED

REMOVABLE INSOLE

LE EXTRAFLEXIBLE PREDISPOSITION
AND SHOCK FOR INSOLE
ABSORBING OUTSOLE

NO INNER

MEDICAL DEVICE

CAL SUITABLE FOR THI CE PREVENTION OF DIABETIC FOOT

ELASTIC

TECNICA 3 KLE BLU

Certification no. 1768378

Color	Blue
Sizes	from 35 to 50

Ref. no. EX012400



SILVER!

Silver Plus silver ion antibacterial treatment, water-resistant and not harmful to the skin.

TECNICA 3 MIK SABBIA

Certification no. 1791287

Color	Sand
Sizes	from 35 to 50

Ref. no. EX012250

TECNICA 3 AGNELLO

Certification no. 1254503

Color	Blue
Sizes	from 35 to 50







TECNICA 3 A KLE NERO

Certification no. 1768377

Color	Black					
Sizes	from 35 to 50					
Ref. no.	EX012600					





Silver Plus silver ion antibacterial treatment, water-resistant and not harmful to the skin.





















REMOVABLE EXTRA FLEXIBLE PREDISPOSITION INSOLE AND SHOCK FOR INSOLE ABSORBING OUTSOLE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

TECNICA 4T

Color	Grey						
Sizes	from 35 to 50						
Ref. no.	EX012050						

TECNICA 4M

Certification no. 983662

Color	Blue
Sizes	from 35 to 50
Ref. no.	EX012050





TECNICA 11

Certification no. 983665

Sizes from 35 to 50	
312es 110111 33 to 30	

Ref. no. EX012200

TECNICA 10

Certification no. 1765775

Color	Grey
Sizes	from 35 to 50
Ref. no.	EX011550





TECNICA T2 AIR

Color	Grey						
Sizes	from 35 to 50						
Ref. no.	EX012350						





















REMOVABLE EXTRA FLEXIBLE PREDISPOSITION INSOLE AND SHOCK FOR INSOLE ABSORBING OUTSOLE

NO INNER STITCHING

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

TECNICA 5T

Certification no. 1511277

Color	Grey
Sizes	from 35 to 50
Ref. no.	EX012350



TECNICA 7

Certification no. 1762241

	Color	Grey
Sizes		from 35 to 50
Ī	Ref. no.	EX012040



TECNICA 6

Color	Grey
Sizes	from 35 to 50
Ref. no.	EX011880























AND SHOCK
ABSORBING OUTSOLE FOR INSOLE

STITCHING

PREVENTION OF DIABETIC FOOT

PLANTARE S-AIR FISIO

Certification no. 1254506

Sizes from 35 to 50 Ref. no. EX010610

PLANTARE S-RIGID FISIO

Certification no. 1254507

Sizes from 35 to 50

Ref. no. EX011150





PLANTARE S-CLOUD FISIO

Certification no. 2142383

Sizes from 35 to 50 Ref. no. EX011490



INTERNAL MEASURES





Lenght








MEASURES OF THE MODELS: 3A KLE NERO, 3 KLE BLU, GOLD 1, GOLD 3, FLEX, 3A, 6, 11, 10, 1, 3T, 3S, 3M, 3E, 3 AIR, 4M, 4T, 3 MIK SABBIA, 3 AGNELLO, SICILIA, AVOLA, SIRACUSA, T2 AIR, 5T, 7, DYNAMIC UOMO, DYNAMIC AIR, SOFT, SASHA, PAUL, GOLD AIR, LARRY, AGRI, ANDRIA.

Sizes	UM	35	36	37	38	39	40	41	42	43	44	45	46	47/48	49/50
Lenght	mm	225	232	239	245	252	259	265	272	279	285	292	299	308	323
Forefoot width	mm	83	84,5	86	87,5	89	90,5	92	93,5	95	96,5	98	99,5	102	105
Heel width	mm	63,5	65	66	67	68	69,5	71	72	73	74	75	77	78,5	81
Forefoot circumference	mm	235	239,5	244	249	253,5	258	262,5	267	271,5	276	280,5	285	292	301
Toe height	mm	35	36	36,5	37	37	37,5	38,5	39	39,5	40	41	41,5	42,5	43,5
Instep circumference	mm	346	357	366	374	383	387	396	403	413	422	431	439	451	470





INTRODUCTION TO THE DIABETIC FOOT BRACES AND POST-OPERATIVE DEVICES LINE

The Tecnica FOOT-ANKLE BRACES and POST-OPE-RATIVE DEVICES are designed for the **treatment of the diabetic foot in its acute phase** and for the pressure unloading in case of **neuropathic and ischemic forefoot or midfoot ulcers**. They are also ideal for **post-surgical treatments** and for all foot surgeries. All uppers have a super-soft padding and no internal

seams in order to minimize the frictions. The multi-velcro closure guarantees a comfortable fit even with bandaged foot and facilitates the adjustment of the footwear. In addition, they have a stabilizing and thermoformable heel-counter and a shock-absorbing rigid and rocking outsole with a non-slip island to ensure better stability.





VARIABLE VOLUME - AMBIDEXTROUS















REMOVABLE INSOLE

FOR INSOLE

NO INNER

DEVICE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

ELASTIC LINING

TECNICA DIAB

Certification no. 983108



Color Black Sizes from S to XXL Ref. no. EX026010



INDICATIONS

- · Diabetic plantar ulcer
- Plantar ulcer
- Post-amputation rehabilitation (up to trans-metatarsal amputation)
- Post-operative rehabilitation of the tarsus
- Fracture of the tarsus
- Management of the chronic phase of Charcot's foot
- Malleolus fracture

CHARACTERISTICS

High leg-foot brace, which is adaptable to the patient, allows the displacement of part of the load on the ankles and calves. The outsole, thanks to the features of movement control and pressure redistribution decompresses the areas affected by the injury, improving thus significantly the ambulation. High thermoformable stiffener. Rigid, anti-flexion and rocking outsole. Variable volume.

Predisposition for the S-Cloud Diab insole. Predisposition for the S-Tal Diab insole.

		Sizes		
S	М	L	XL	XXL
34-35-36	37-38-39	40-41-42	43-44-45	46-47-48

TECNICA POST

Certification no. 983646

Color Black Sizes from S to XXL Ref. no. EX024580





INDICATIONS

- Diabetic plantar ulcer with vascular complications affecting the leg (combined with the S-Cloud Diab insole)
- · Neuroischemic ulcer
- Management of the injured foot or foot with ulcerative risk
- Compound fracture of the metatarsals
- Post-amputation rehabilitation of the toes
- Fracture of the toes
- Surgical treatment of the toes
- Pre-surgery stabilization

CHARACTERISTICS

Low leg-foot brace, which is adaptable to the patient. The outsole, thanks to the features of movement control and pressure redistribution decompresses the areas affected by the injury, improving thus significantly the ambulation and the recovery time.

High thermoformable stiffener.

Rigid, anti-flexion and rocking outsole. Variable

Predisposition for the S-Cloud Diab insole. Predisposition for the S-Tal Diab insole.

		Sizes		
S	М	L	XL	XXL
34-35-36	37-38-39	40-41-42	43-44-45	46-47-48







VARIABLE VOLUME - AMBIDEXTROUS

















REMOVABLE **INSOLE**

FOR INSOLE

NO INNER

MEDICAL

SUITABLE FOR THE PREVENTION OF **DIABETIC FOOT**

UPPER

PADDED LINING

TECNICA TERAP



TECNICA HEEL Certification no. 1294420 Color Black Sizes from S to XL



INDICATIONS

- · Diabetic ulcerations
- Plantar ulcerations
- Neuroischemic ulcerations
- Fracture of the metatarsals
- Surgery and fractures of the toes
- · Rehabilitation after toes amputation
- Pre-surgery stabilization
- Instability of the hindfoot

CHARACTERISTICS

Easy fit and maximum safety for the patient during the acute phase. The thermoformable upper in pel-flex, which is partially waterproof, guarantees thermal and hygienic protection to the foot. Radiolucent and anti-friction. Thermoformable stiffener.

Rigid, anti-flexion and rocking outsole. Variable volume. Predisposition for the insertion of the S-Cloud Diab insole.

		Sizes		
S	М	L	XL	XXL
34-35-36	37-38-39	40-41-42	43-44-45	46-47-48

INDICATIONS

- Post-surgery rehabilitation of the heel area
- · Wounds or ulcers of the heel area
- Post-surgery fractures and/or calcaneus ulcers (second treatment) and ulcers at the insertion of the Achilles tendon
- Bursitis

CHARACTERISTICS

The balance of the weight is moved to the front of the foot thanks to the unloading insole and thanks to the outsole, which is designed to totally release the heel from the pressure load. The square shape of the outsole offers more space for bandages. Movable rear-cover to preserve the heel from impacts. Non-slip outsole to ensure maximum safety during walking. Total opening with velcro-straps. Removable insole.

	Siz	es	
S	М	L	XL
34-35-36	37-38-39	40-41-42	43-44-45

TECNICA CT DIAB

Certification no. 1255252

Color	Black
Sizes	from S to XXL
Dof no	EVANADO



PLANTARE S-CLOUD DIAB

from S to XXL

Certification no. 1903610

Sizes

Ref.	EX020930	
no.		
	000	
teci	nica)	

PLANTARE S-TAL DIAB

Certification no. 1255254

Sizes	from S to XXL
Ref.	EX020970







VARIABLE VOLUME - AMBIDEXTROUS

















REMOVABLE INSOLE

PREDISPOSITION FOR INSOLE

NO INNER

MEDICAL DEVICE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

ELASTIC UPPER

PADDED LINING

ROCKING OUTSOLE

TECNICA JOLLY

Certification no. 983648



Color Black
Sizes from S to XXL

Ref. no. EX023030



TECNICA CT-DIAB

INDICATIONS

- •Post-surgery of hallux valgus
- Hammer toes
- Morton's neuroma
- Fracture of the toes
- Nail removal
- Plantar wart removal

CHARACTERISTICS

Rehabilitative foot orthosis, with radiolucent upper, which is antifriction and washable. The non-slip, rigid and rocking outsole unloads the pressure from the forefoot, giving thus greater stability without altering the posture, ensuring a rapid recovery of all functions.

Thermoformable stiffener. Variable volume. Predisposition for the S-Cloud Diab insole. Predisposition for the S-Tal Diab insole.

		Sizes		
S	М	L	XL	XXL
34-35-36	37-38-39	40-41-42	43-44-45	46-47-48

TECNICA VALGUS

Certification no. 983649

Color Black
Sizes from S to XXL
Ref. no. EXO22540



INDICATIONS


- Post-surgery of hallux valgus, hammer toes, Neuroma of Morton
- Fracture of the toes, nail removal
- Plantar wart removal
- Best solution for diabetic feet and indicated in presence of trans-metatarsal amputations
- Ray amputations
- Chronic Charcot of the forefoot and of the midfoot
- Minor amputations of the forefoot


CHARACTERISTICS

The rigid and rocking outsole cushions the step and immobilizes the metatarsophalangeal area, reducing the pressure, allowing thus a total discharge of the trans-metatarsal area. Allows to accommodate deformities and allows to avoid dangerous frictions. The wide plant and the anti-stumbling outsole ensure maximum stability during walking. Higher stabilizing stiffener for better foot support. Variable volume. Removable insole.

		Sizes		
S	М	L	XL	XXL
34-35-36	37-38-39	40-41-42	43-44-45	46-47-48

INTERNAL MEASURES









MEASURES OF THE MODELS: DIAB, POST, TERAP, JOLLY, VALGUS.

		, - ,	,	,		
Sizes	UM	S	М	L	XL	XXL
Lenght	mm	230	250	270	290	310
Forefoot width	mm	95	100	105	110	115
Heel width	mm	78	80	82	86	90





VARIABLE VOLUME - AMBIDEXTROUS

















REMOVABLE INSOLE

PREDISPOSITION FOR INSOLE

NO INNER

MEDICAL DEVICE

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

ELASTIC UPPER

PADDED LINING

UTSOLE

TECNICA STELT

Certification no. 2142377

Color Black/White
Sizes from 35 to 48
Ref. no. EX022980



TECNICA CT STELT

Certification no. 2142382

Color	Black
Sizes	from 35 to 48
Ref. no.	EX020930



PLANTARE S-TAL STELT

Certification no. 2142380

Sizes	from 35 to 48
Ref.	EX020970
no.	



INDICATIONS

- · Rigid biomechanical outsole
- 20mm outsole which minimizes discrepancy with the contralateral limb
- S-Tal Stelt insole included
- Stabilizing and thermoformable stiffener
- Multi-velcro closure
- Ambidextrous
- Predisposition for the insertion of the S-Cloud Stelt unloading insole

CHARACTERISTICS

Tecnica STELT is an ambidextrous footwear, indicated both as a brace and as a post-operative device. It is particularly suitable for the relieve of pressure in case of neuropathic and ischemic ulcers in the forefoot and midfoot. Thanks to the rigid outsole it guarantees the maximum stability and immobilization of the foot, optimizing thus the recovery times. It is also used as a postoperative device for all surgeries on the forefoot. Thanks to its very wide fit, the footwear can be used also with important bandages. The height of the outsole of 20 mm allows to maintain the physiological posture and to minimize the discrepancy with the contralateral limb, helping thus to simulate a natural walk. The S-Tal Stelt insole is included in the package and has an 8° talism angle, which allows to transfer the pressure from the forefoot to the hindfoot.

			Sizes			
35-36	37-38	39-40	41-42	43-44	45-46	47-48



PLANTARE S-CLOUD STELT

Certification no. 2142384

Sizes	s from 35 to 48
Ref.	EX020930
no.	YXX
	tecnica





VARIABLE VOLUME - AMBIDEXTROUS



STITCHING









MEDICAL SUITABL DEVICE PREVEI

SUITABLE FOR THE PREVENTION OF DIABETIC FOOT

PADDED LINING

ROCKING OUTSOLE

TECNICA BIOMED

Certification no. 1765774

Color Black/White
Sizes from 35 to 46

Ref. no. EX022100



PLANTARE S-CLOUD BIOMED





		Sizes POS	T-OPERATIVE		
35-36	37-38	39-40	41-42	43-44	45-46



TECNICA 20 RIGID

Certification no. 983651

Color	Black
Sizes	from 35 to 46
Ref. no.	EX021930

INDICATIONS

• II° post-surgical treatment

CHARACTERISTICS

Super light and very stable footwear thanks to the thermoformable stiffener, with open front, particularly suitable for deformities and injuries in the forefoot area and ideal also for physiotherapy for functional recovery. The variable internal volume of the upper, which can be completely opened, facilitates the fit of bandaged or painful feet in the phase following the brace. The stiffener is customizable by the orthopedic technician so that the fit can be adapted to the patient's foot. Non-slip outsole for maximum safety during walking. Breathable upper in washable microfiber.



		Sizes POST	T-OPERATIVE		
35-36	37-38	39-40	41-42	43-44	45-46

INTERNAL MEASURES









MEASURES OF THE MODELS: BIOMED, STELT, 8, 16, 20 RIGID.

Sizes	UM	35/36	37/38	39/40	41/42	43/44	45/46	47/48
Lenght	mm	235	250	260	270	285	295	310
Forefoot width	mm	85	90	95	100	105	110	115
Heel width	mm	72	74	74	78	80	82	84





VARIABLE VOLUME - AMBIDEXTROUS









NO INNER

SUITABLE FOR THE PREVENTION OF

LINING

TECNICA 16

Certification no. 983652

Color	Grey
Sizes	from 35 to 48
Ref no	EY021350



INDICATIONS

- Post-surgical treatment of forefoot fractures
- Hallux valgus, hammer toes
- Bursitis
- Forefoot and interdigital ulcerations

CHARACTERISTICS

Designed in order to provide protection and relief to the forefoot. The outsole in talism release the operated area from the pressure which is shifted to the metatarsal and to the heel. Allows a precocious walking, anticipating functional recovery. Wide shaped outsole and rigid stiffener for maximum stability. Variable volume. Upper in breathable microperforated fabric washable at 40°.

		Size	s POST-OPE	RATIVE		
35-36	37-38	39-40	41-42	43-44	45-46	47-48

TECNICA 8

Certification no. 983650

Color	Grey
Sizes	from 35 to 48
D-f	EV/0040E0



TECNICA MOVIDA

Certification no. 1255251

Color	Blue
Sizes	from 35 to 46
- C	=>/0000=0



INDICATIONS

- Post-surgical treatment of forefoot fractures
- · Hallux valgus, hammer toes
- Bursitis
- Forefoot and interdigital ulcerations

CHARACTERISTICS

Designed in order to provide protection and relief to the forefoot. The outsole in talism release the operated area from the pressure which is shifted to the metatarsal and to the heel. Allows a precocious walking, anticipating functional recovery. Wide shaped outsole and rigid stiffener for max mum stability. Variable volume. Upper in breathable microperforated fabric w shable at 40°. With front rim in order to apply the Tecnica C8 toe-cover for thermal and hygienic protection.

TECNICA C8

Certification no. 983666 Ref. no. EX020440



		Sizes	POST-OPE	RATIVE		
35-36*	37-38*	39-40*	41-42**	43-44**	45-46**	47-48**

*	Size	Small
**	Si70	Largo

INDICATIONS

• Post-surgical treatment

CHARACTERISTICS

Cast shoe suitable for cast bandages. Upper in breathable jeans fabric with laces and robust non-slip outsole, with a 20mm high peripheral edge for protection from water.

Sizes TECNICA MOVIDA				
35-36-37	38-39-40	41-42-43	44-45-46	



















THE FIRST SAFETY & ORTHOPAEDIC FOOTWEAR





THE PROTECTION FOR YOUR SENSITIVE FEET

























INTRODUCTION

EXTREME LINE

Footwear with double certification

From the technology and experience of ExtroStyle, a leading manufacturer of footwear for foot health, comes a unique line of footwear: **EXTREME!** Especially designed to protect sensitive feet during work. **THIS FOOTWEAR POSSESSES A DOUBLE CERTIFICATION: SAFETY CERTIFICATION AS IT IS A PPE PERSONAL PROTECTIVE DEVICE ACCORDING TO EN ISO**

20345:2011 AND, AT THE SAME TIME, AN ORTHOPAEDIC FOOTWEAR MD MEDICAL DEVICE ACCORDING TO CE MDR 745/2017. They have been created for women and men who are active every day in the world of work wearing a protective device, with the opportunity to wear a shoe designed for the insertion of a customised foot orthosis.

PPE - PERSONAL PROTECTIVE EQUIPMENT

















ANTI-PUNCTURE INSERT

HEEL AREA ENERGY ABSORPTION

ANTISTATIC

DYNAMIC WATERPROOFING OF THE UPPER

RESISTANCE TO OILS AND HYDROCARBON

SLIPPING RESISTANCE

200-JOULE RESISTANT TOE CAP

UPPER

EXTREME GREY (flex outsole)

Sizes available from 36 to 48 Certification no. 2207738 REF. N° EX067250







EXTREME WHITE (flex outsole)

Sizes available from 36 to 48 Certification no. 2207903 REF. N° EX067250







The EXTREME outsole has a special tread design that guarantees high anti-slip properties, both in the longitudinal and transverse direction, as well as being antistatic and resistant to hydrocarbons, oils and greases in general.



















DESCRIPTION

EXTREME LINE AS PERSONAL PROTECTIVE EQUIPMENT

EXTREME safety footwear for professional use is certified in accordance with the latest regulations; it is suitable for most types of use as it is S3 SRC PROTECTION grade. It is a soft and a lightweight shoe, with leather upper and technical fabric inserts, anti-crush toe cap, antistatic and multilayer ballistic fabric insert with the absolute guarantee of a puncture barrier, which is also highly flexible.

TEPOR's waterproof and breathable membrane does not allow water to enter the shoe and, at the

same time, allows the foot's moisture to escape. **TE-POR** thereby guarantees maximum comfort while keeping the foot dry with an optimal microclimate inside the shoe.

The **polyether sole** guarantees absolute comfort due to its lightness, ensures flexibility and absorption of vibrations and shocks, is antistatic and resistant to hydrocarbons and oils and greases in general. The extreme compactness and the special tread design guarantee excellent anti-slip properties in both longitudinal and transverse directions.

EN ISO 20345:2011 STANDARD SAFETY CATEGORY				
MINIMUM REQUIREMENTS	REFERENCE	S 3		
200-JOULE RESISTANT TOE CAP	S	A STATE OF THE STA		
RESISTANCE TO OILS AND HYDROCARBONS	FO	Cuthan.		
ENERGY ABSORPTION IN THE HEEL AREA	E	O REPUBLICATION OF THE PROPERTY OF THE PROPERT		
ANTISTATIC	A	- International Control		
ANTI-PUNCTURE	Р	- International Control		
DYNAMIC WATERPROOFING OF THE UPPER	WRU			
SLIPPING RESISTANCE	SR	Onton,		

INDICATIONS FOR USE AS PPE

The S3 standard corresponds to work shoes with full protection, which is why they are in high demand. They combine the properties of S2 with the anti-puncture layer typical of S1 footwear. An ideal combination for many high-risk work environments, such as:

LIGHT INDUSTRY • Food and Beverage • Pharmaceutical • Mechanical • Chemical • Automotive • Textile • Footwear • Clothing • Wood • Paper • Construction.

HEAVY INDUSTRY •• Oil and Gas • Metallurgy • Mechanical carpentry • Shipyards.

TRADES • Personal Care • Electrician • Carpenter • Mechanic • Plumber • Mason.

LOGISTICS AND TRANSPORT • Warehousing • Agriculture and livestock • Security and Surveillance.

CATERING AND HOSPITALITY • Hospital Healthcare • Services and Tertiary Sector.







INDICATIONS FOR

REMOVABLE INSOLE

The standard removable insoles used for the **EXTREME** line are composed of a high-density, open-cell foam that provides total breathability while protecting against perspiration due to its moisture-absorbing and moisture-wicking qualities, keeping the foot ventilated, dry and odour-free, and providing a pleasant sensation during walking. The ergonomic design of the insole provides a high degree of comfort, minimising strain and tension in the plantar support, and cushioning impact pressure both while walking and also in a static position.

It is guaranteed that the performance has been determined by carrying out technical tests on the footwear fitted with a **removable insole**.

If it is necessary to replace the standard removable insole included with the footwear at the time of purchase, it must be replaced with an identical one supplied by the manufacturer in order to not alter the certified configuration. If it is necessary to insert a corrective and/or antalgic foot orthosis for preventive or curative foot care, in this case, the treatment shall only be performed by a qualified orthopaedic technician, who must comply with the manufacturer's recommendation that the maximum height of the orthosis near the toe cap does not exceed 6 (six) millimetres.

Compliance with the manufacturer's instructions is indispensable in order to not alter the technical tests carried out during certification.



















DESCRIPTION

EXTREME S3 SRC





The EXTREME heel support is padded and robust, it gently embrace the heel area while providing stability and protection. The 1.6-1.8 cm thick padded leather upper, with almost no seams, accommodates the foot in a protective and shockabsorbing envelope.

The highly breathable AIRTEK **lining** provides the inside of the shoe with a complete air exchange with a "chimney effect"



cork for lightness and thermal isolation.

The ergonomic design of the insole provides a high degree of comfort, minimising strain and tension in the plantar support, and cushioning impact pressure both while walking and also in a static position.

flexible.

THE EXTREME LINE IS METAL FREE

The absence of metal parts does not affect the safety and the protection of the wearer; on the contrary, it guarantees very high safety standards during any kind of performance. EXTREME Metal Free shoes are appreciated by all those who work in environments with the presence of metal detectors, such as: airports, banks, public offices, post offices, museums, boarding areas, customs offices and much more.





















DESCRIPTION

EXTREMELINEAS A MEDICAL DEVICE

Medical Device

The **EXTREME** line is especially designed to protect sensitive feet during work activities. It is produced in two versions with a FLEXIBLE outsole and with a RIGID outsole in accordance with international protocols for the treatment and care of diabetic feet. The 1.6–1.8 mm thick full-grain leather upper combined with the waterproof and breathable **TE-POR** membrane and the highly breathable **AIR-TEK** lining provide the inside of the shoe with a complete "chimney effect" air exchange while maintaining an

optimal microclimate. The absence of internal seams ensures that there are no chafing and pressure points. It envelops the anatomy of the foot and any deformities while providing comfort and preventing annoying friction. The biomechanical **EXTREME** outsole allows a natural transition between the three phases of walking: it ensures a reduction of pressure peaks on the forefoot, limits wear and tear on the joints, enables spontaneous walking and reduces the sensation of fatigue.

MD - MEDICAL DEVICE



REMOVABLE INSOLE



MAXIMUM SHOCK-ABSORBING CAPACITY EXTRA-FLEX OUTSOLE



PREDISPOSITION FOR INSOLE



NO INTERNAL SEAMS



MEDICAL DEVICE



SUITABLE FOR THE PREVENTION OF DIABETIC FOOT DISEASE



PADDED LIPPER





EXTROSTYLE S.R.L. IN ITS ROLE AS MANUFACTURER, DECLARES THAT THE EXTREME BRAND FOOTWEAR LINES ARE CLASS I NON-INVASIVE MEDICAL DEVICES, COMPLYING WITH THE REQUIREMENTS OF EU MDR 745/2017.



















DIRECTIONS FOR USE

EXTREME FLEXIBLE OUTSOLE

EXTREME FLEXIBLE OUTSOLE with removable insole is a medical device used for restraint purposes with a customised insole made by qualified health professionals for the treatment of alterations in the morphology of the foot. Suitable for all those who use a customised foot orthosis to correct, prevent or treat various temporary or chronic foot disorders. Particularly indicated for the treatment of rheumatic

foot and grade 0 and 1 prevention of diabetic foot disease also in presence of neuropathies accompanied by vasculopathies, as well as metatarsalgia, polyarthritis and arthrosis. The biomechanical sole facilitates the step rolling by supporting the foot without constricting or forcing it. It reduces pressure peaks on the forefoot and limits the load on the joints.

INDICATIONS FOR USE IN THE PRESENCE OF:

- Primary prevention in diabetic foot disease with risk grade 0 or 1
- Polyarthritis
- → Morton's neuroma
- → Arthrosis
- Metatarsalgie


EXTREME RIGID OUTSOLE

EXTREME RIGID OUTSOLE is a medical device used for secondary prevention for diabetic patients with risk levels 2 and 3, where there is the need to relieve plantar pressure: designed for those suffering from previous ulcers or amputations, it is also suitable for patients suffering from rheumatism

and severe deformities. The rigid, biomechanical outsole limits joint stresses of the foot during walking, by significantly reducing pressure peaks in the forefoot and redistributing them correctly, thus reducing the risk of reulceration in diabetic patients by more than 60%.

INDICATIONS FOR USE IN THE PRESENCE OF:

- Diabetic patient with risk level 2
- Previous ulcers at risk of recurrence with risk grade 2
- Deformity, ulceration or previous amputation with risk grade 3
- Finger and tarsal arthrodesis
- → Charcot foot





















EXTREME LINE

INTERNAL MEASURES

FIT 12 UM mm						
SIZE	INSTEP CIRCUMFERENCE	LENGTH	FOREFOOT CIRCUMFERENCE	FOREFOOT WIDTH	HEEL WIDTH	TOE HEIGHT
48	294,5	325	296	101	71	36,5
47	290	318	291	99,5	70	36
46	285,5	311	286	98	69	35,5
45	281	304	281	96,5	68	35
44	276,5	298	277	95,5	67	34,5
43	272	291	272	94	66	34
42	267,5	284	268	93	65	32,5
41	258	278	258	89	62	32
40	293,5	271	253	87,5	61	31,5
39	249	264	249,5	86	60	31
38	244,5	258	244	84	59	30,5
37	240	251	239,5	83	57,5	30
36	235,5	244	235,5	81,5	57	29,5
35	231	238	231	80	56	29
Forefoot width Heel width Forefoot circumference Forefoot circumference Toe height						

















NOTES		
<u></u>	 	 ······
	 	 ·····

DATA TREATMENT

We inform you that, through the fulfillment of your order, it will be necessary to process the data you transmitted to our structure in electronic, magnetic, and paper format. The ExtroStyle company srl declares to be in compliance with the security measures provided by current legislation in the field of Privacy, Legislative Decree number 196 of the 30th June 2003 and the supplementary rules provided by the G.D.P.R (EU regulation 679/2016).





TOOLS FOR THE STORE

EXTRO STYLE PROVIDES STORES WITH A WIDE CHOICE OF DISPLAY AND SUPPORT MATERIAL







FOOT MEASURING DEVICE

STORE DISPLAY



COUNTER DISPLAYS





CATALOG

THEMATIC CATALOGS



SHOE STAND



NOTES







NOTES

		······································





NOTES









